
# WILOVERT

INVERSE MIKROSKOP-SERIE

INVERTED MICROSCOPES

## INVERSE MIKROSKOPE

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## INVERTED MICROSCOPES

For a wide variety of applications, inverted microscopes are an indispensable tool for the investigation of living material. On the other hand, they have become highly specialized systems for certain routine applications and research fields. This is exactly what the Wilovert line of microscopes is designed for. Its components include:

### MICROSCOPE STAND

The U-shaped stand of the Wilovert series is built from one single cast part as an optimal prerequisite for high stability. Micromanipulators can be adapted to the Wilovert microscopes with ease.

### LIGHT SOURCES

For transmitted-light brightfield applications, the Wilovert carries a 12 V/30 W halogen lamp. For incident-light fluorescence investigations, an HBO 100 W mercury high-pressure lamp can be adapted.

### OBJECTIVES

As the specimens are contained in culture vessels of different materials, heights and thicknesses at bottom, inverted microscopes require dedicated objectives. The Helmut Hund GmbH offers objective/condenser combinations with long free working distances.

### CONDENSERS

The available condensers are:

- Condenser NA 0.25 for transmitted-light brightfield with a free working distance of 58 mm
- Combination condenser NA 0.25 for transmitted-light brightfield and phase contrast with a free working distance of 58 mm
- Combination condenser NA 0.4 for transmitted-light brightfield and phase contrast with a free working distance of 38 mm

### HOLDERS FOR CULTURE VESSELS

Culture vessels like Petri dishes, tissue culture vessels, microtest trays and microtiter trays are inserted in special holding frames mounted to the stage of the Wilovert microscopes. The object guide ensures that specimens can be scanned completely.

**The configurations contained in this brochure are only examples for the possible variants of the Wilovert.**

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### WILOVERT STANDARD HF

Standard configuration for routine investigations of transparent objects with high-contrast structures that can be observed in brightfield. The advantage of this model is its simplicity and ease-of-use that permits ample room for your applications.

- Eyepieces WF 10x/18
- Achromatic objectives  
A 4/0.10  
A 10/0.25

### WILOVERT STANDARD PH 20

Routine microscope for the observation of cell cultures in brightfield and phase contrast.

- Eyepieces WF 10x/18
- Achromatic and semiplan-achromatic objectives  
A 4/0.10  
A Ph 10/0.25  
SPL Ph 20/0.35 LD
- Combination condenser NA 0.25

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### WILOVERT STANDARD HF 40

Standard configuration for the observation of cell cultures in transmitted-light brightfield.

- Eyepieces WF 10x/18
- Achromatic and semiplan-achromatic objectives  
A 4/0.10  
A Ph 10/0.25  
SPL 20/0.35 LD  
SPL 40/0.60 LD in correction mount
- Condenser NA 0.25
- Object guide with left-hand drive

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### WILOVERT STANDARD PH 40

Standard configuration for the observation of cell cultures in transmitted-light brightfield and phase contrast.

- Eyepieces WF 10x/18
- Achromatic and semiplan-achromatic objectives  
A 4/0.10  
A Ph 10/0.25  
SPL Ph 20/0.35 LD  
SPL Ph 40/0.60 LD in correction mount
- Combination condenser NA 0.25
- Object guide with left-hand drive

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## WILOVERT AFL FOR INCIDENT-LIGHT FLUORESCENCE

Fluorescence plays an important role in the development of investigation and measurement techniques in molecular biology and biochemistry. For high-contrast and brilliant fluorescence images, the proven HBO 100 W mercury high-pressure lamp is employed with the Wilovert AFL. Filter slides for two or four filter sets allow rapid switching between investigations with different fluorochromes. The Hund filter programme contains a long line of sets for standard applications.

### WILOVERT AFL 40

For investigations in transmitted-light brightfield and incident-light fluorescence (e. g. for HLA typing).

- Eyepieces WF10x/18
- Achromatic and semiplan-achromatic objectives  
A 4/0.10  
A 10/0.25  
SPL 20/0.35 LD  
SPL 40/0.60 LD in correction mount
- Condenser NA 0.25
- Filter set for FITC
- HBO 100 W lamp
- Object guide with left-hand drive

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### WILOVERT AFL 40 PHAKO

Universal configuration for investigations in incident-light fluorescence in combination with transmitted-light phase contrast.

- Eyepieces WF 10x/18
- Achromatic and semiplan-achromatic objectives  
A 4/0.10  
A Ph1 10/0.25  
SPL 20/0.35 LD  
SPL Ph2 20/0.35 LD  
SPL Ph2 40/0.60 LD in correction mount
- Combination condenser NA 0.25
- Filter set for FITC
- HBO 100 W lamp
- Object guide with left-hand drive

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### OBJECTIVES

The objective programme includes:

- Achromatic objectives
- Semiplan-achromatic objectives
- Phase-contrast achromatic objectives
- Phase-contrast semiplan-achromatic objectives

All objectives are designed for a tube length of 160 mm.

### EYEPIECES

Wide-field eyepieces with magnifications of 10x and 12.5x are standard for the Wilovert. Eyepieces for spectacle-wearers are also available. For microscopic measurements, we recommend measuring and counting eyepieces with graticules.

### TRINOCULAR PHOTOTUBES

#### Trinocular Tube 70/30

70 % of the light are directed to the camera port, 30 % are available for a simultaneous visual observation.

#### Trinocular Tube 100/100

100 % of the light are directed either to the camera port or to the eyepieces. This system is most useful for all applications where only small light intensities are available, e. g. for fluorescence microscopy.

Cameras are fitted to the camera port via C-mount adaptors or dedicated digital camera adaptor sets.

### MICROMANIPULATORS AND MICROINJECTORS

For working with living cells, e. g. for the injection of substances into the cell nucleus or into the cytoplasm, micro-manipulators and microinjectors are available. For these applications, the microscope stage should be vibration-free to allow microinjections and electrophysiological investigations. To this end, the stage is screwed to both sides of the microscope stand. As focusing is not done by adjusting the stage height, but by adjusting the objective revolver, the stage-specimen-microinjector system remains permanently stable.

# hund

Helmut Hund GmbH  
Wilhelm-Will-Str. 7  
D-35580 Wetzlar  
Telefon: +49 (0) 6441 2004-0  
Telefax: +49 (0) 6441 2004-44  
E-Mail: [zentrale@hund.de](mailto:zentrale@hund.de)  
Internet: [www.hund.de](http://www.hund.de)