

																	
	medicus plus	medicus plus PH	H 600 Wilo-Prax Achro	H 600 Wilo-Prax PL	H 600 Wilo-Prax ICS	H 600 Wilozyt Plan	H 600 PH Plan	H 600 LED AFL Myko	H 600 HP LED	H 600 LL HP 100	H 600 AFL Achro	H 600 AFL Plan	H 600 AFL PH 40	Wilovert Standard HF	Wilovert Standard PH	Wilovert AFL	
Hematology	●	●	●	●	●	●											
Urology		●				●	●										
Gynecology						●	●								●		
Dermatology							●	●							●		
Pathology, Bacteriology				●	●		●							●	●	●	
Immunology								●			●	●	●				●
Detection of Parasites								●	●	●				●	●	●	

	Med-Prax 3	medicus plus	medicus plus PH	H 600 Wilo-Prax Achro	H 600 Wilo-Prax PL	H 600 Wilo-Prax ICS	H 600 Wilozyt Plan	H 600 PH Plan	H 600 Wilo-Brau	H 600 BS	H 600 LED AFL Myko	H 600 HP LED	H 600 LL HP 100	Wilovert Standard HF 20	Wilovert Standard HF 40	Wilovert Standard PH 20	Wilovert Standard PH 40	Wilovert AFL 20	Wilovert AFL 40
	Stained Tissue Sections	●	●		●	●	●												
Fungi and Mycoses			●				●	●			●								
Breweries									●										
Sewage Plants										●									
Plankton, Microbes														●	●				
Sediments														●	●				
Cell Cultures														●	●	●	●	●	●
Tissue Cultures														●	●				
Live Blood												●	●						

	Med-Prax 3	medicus plus (PH)	H 600 AM 50	H 600 AM AL/DL 50	H 600 POL, LL POL 50	Wiloskop ST	Wiloskop ST-AD	Wiloskop ST-S
Metallurgy			●	●	●			
Material Testing			●	●	●	●	●	●
Quality Assurance, MFG			●	●	●	●	●	●
Section						●	●	●
Education	●	●						

Remarks:

- All H 600 equipments available with Plan ICS optics (infinity corrected) available upon request.
- Simple and sophisticated polarization units available for all upright microscopes.
- H 600 AFL equipped with 50-W HBO illumination, Wilovert AFL with 100-W HBO illumination.
- All microscopes (except for Med-Prax 3) can be equipped with trinocular tubes upon request. H 600 HP LED and H 600 LL HP 100 are already equipped with trinocular tubes.